

# Separations | HPLC Columns ProteCol™ Amino

## Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

### Features and Benefits

Allows basic compound separation in normal phase and carbohydrate analysis.

### Recommended Applications

Reverse phase carbohydrate analysis and normal phase basic compound separations.

### Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	USP Code
Amino	120	5	350	4.0	L8

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).

