

Separations | HPLC Columns ProteCol™ Cyano

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and Benefits

Provides chromatographic retention both in normal and reversed phase separation due to its moderate polarity.

Recommended Applications

Separation of flavonoids, extraction of polar compounds from non-polar samples as well as the analysis of samples containing analytes with a wide range of hydrophobicity.

Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	USP Code
Cyano	120	5	300	5	L10

For more information about this product visit www.sge.com or contact techsupport@sge.com.

