

# Separations | HPLC Columns ProteCol™ Phenyl Hexyl

## Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

### Features and Benefits

ProteCol Phenyl-Hexyl uses a hexyl-linked phenyl phase where the hexyl alkyl chain delivers unique selectivity and increased hydrolytic stability when compared to propyl-linked chemistry.

### Recommended Applications

Recommended for compounds that are difficult to separate on standard phenyl phase columns.

### Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	pH	USP Code
Phenyl Hexyl	120	3, 5	350	10.0	1-10	L11

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).

