

# Separations | HPLC Columns ProteCol™ SCX

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

## Features and Benefits

- High loading capacity and pressure limit
- Robust bonding technology
- High density bonding
- Bonding aromatic sulfonic acid groups

## Recommended Applications

Small organic bases.

## Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

PEEK coated column hardware available to reduce non-specific interactions when analyzing compounds with metal chelating activity.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	pH
SCX	120	3, 5	350	1-10

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).

