

Separations | HPLC Columns ProteCol™ Ultra

Turn Imagination into Results

Complementing the ProteCol™ analytical range. These columns use 1.8 µm particles for UHPLC use.

Features and Benefits

Specifically designed for UHPLC use, sub 2 micron phases enable use up to 19000 psi.

Recommended Applications

UHPLC.

Product Specifications

The following phases are available in 2 mm ID columns, with 50, 100, and 150 mm lengths:

- ProteCol Ultra C18
- ProteCol Ultra C8
- ProteCol Ultra Amino
- ProteCol Ultra Cyano
- ProteCol Ultra HILIC FL
- ProteCol Ultra HILIC PI
- ProteCol Ultra PFP
- ProteCol Ultra Phenyl
- ProteCol Ultra Phenyl Hexyl
- ProteCol Ultra Polar
- ProteCol Ultra SCX

For more information about this product visit www.sge.com or contact techsupport@sge.com.

